

SAP / VID Bracket Placement

	UR7	UR6	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5	UL6	UL7	DIFF
SAP	3.0	3.5	4.0	4.5	5.5	6.0	6.5	6.5	6.0	5.5	4.5	4.0	3.5	3.0	3.5
SAP +	3.0	3.5	4.0	5.0	6.0	6.3	7.0	7.0	6.3	6.0	5.0	4.0	3.5	3.0	4.0
SAP ++	3.0	3.5	4.5	5.5	6.5	6.5	7.5	7.5	6.5	6.5	5.5	4.5	3.5	3.0	4.5

	LR7	LR6	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5	LL6	LL7	DIFF
SAP	5.0	5.0	5.0	5.5	5.0	4.0	4.0	4.0	4.0	5.0	5.5	5.0	5.0	5.0	-1.0
Deep Bite	5.0	5.0	5.0	5.0	4.5	3.5	3.5	3.5	3.5	4.5	5.0	5.0	5.0	5.0	-1.5
Open Bite	5.0	5.0	5.0	5.5	5.5	5.0	5.0	5.0	5.0	5.5	5.5	5.0	5.0	5.0	0

Bracket Gauges

Provided Courtesy of Dr. Tom Pitts

Measured to Bracket Slot

Copyright © 2019, RGB Ortho Systems LLC, All Rights Reserved

